


# E-Wave

## NeuroMuscular Electrical Stimulator



Rx Only - Made in America

### Adjustable Parameters

Amplitude:	0-100 mA
Frequency:	4-100 Hz
Pulse-Width:	50-300 $\mu$ sec
Ramp Up:	0.1 to 60 Sec.
On-Time:	0.5 to 30 Sec.
Ramp Down:	0.1 to 60 Sec.
Off-Time:	0.1 to 60 Sec.
Waveforms:	Symmetrical Biphasic (AC) or Monophasic (DC)
TX Timer:	Continuous or 10-100 Minutes

### Compatible Accessories

The remote switch allows either the patient or therapist the ability to trigger a muscle contraction manually.



External Hand Switch



External Heel Switch

### Indications For Use

- Prevention or retardation of disuse atrophy
- Muscle re-education
- Relaxation of muscle spasms
- Increase local blood circulation
- Maintain or increase range of motion

