



Our Innovative design makes it possible for patients to go home with the same quality traction they receive during physical therapy.



SPECIFICATIONS:

Maximum force: 50 lbs.

Angles of incline: 10, 15, and 20 degrees


Hand pump: Displays force in both pounds and kilograms

Neck wedges: independently adjustable

Headrest material: Memory foam

Testing: 15,000 cycles which equals using 3 times daily for 15 years

Fully Adjustable Angle of Incline without the need to order extra parts




Hand Pump
our patented hand pump makes it easy to increase or decrease traction force

Quick Release Button
push to release the traction force

Memory Pillow
conforms comfortably to the head

Independently Adjustable Neck Wedges
provide a custom fit



More Information:

Seneca Sadler
1-800-495-6670 x2333
ssadler@zynex.com

Available to your patients for home use with a prescription via Zynex Medical's EZ-Prescribe Program

Distributed by:

