

Zynex Pro Post-op Knee Brace

Protection & Stability

The Zynex Pro Post-Op Knee Brace is designed to provide protected, controlled range-of-motion (ROM) for patients recovering from knee surgery or those who have knee injuries or instabilities.

Streamlined Design

- Telescoping struts for length adjustment from 17"- 29" (43.18 cm - 76.33 cm)
- Simple rotate and release Flexion/Extension settings
- Unique safety lock to increase patient compliance
- 100% tool-free adjustments
- Quickly locks at 0°, 15°, 30°, 45°



Sizing Options

Standard - Up to 29" (76.33 cm)
X-Large - Up to 35" (88.9 cm)

Indications for Use

- Post-op rehabilitation
- Soft tissue knee injuries
- Knee sprains/strains
- Tibial plateau fractures
- Osteochondral repairs
- Meniscal repairs
- Patella tendon repairs
- Condylar fractures
- ACL, PCL, MCL, LCL repairs/injuries
- HTOs (High Tibial Osteotomy)
- Stable or internally fixed fractures of the knee

